



PARTS CATALOGUE





Cables and Cylinders







6540 Tomken Road
Mississauga ON L5T 2E9
905.564.2250
905.564.5424
1.800.362.0873
www.mapleleafwheelchair.ca





Cables and Cylinders







Part No.	Part Name	Part Description	Cost	
TRP004 G200	200 N cylinder Long Sold with [CC009]. (SUSPA cylinder Trigger)	Long Suspa cylinder having a 200 N load. -Used in both Refresh Commode Cylinders.	240/ each	
TRP004 G400	400 N cylinder Long Sold with [CC009]. (SUSPA cylinder Trigger)	Long Suspa cylinder having a 400 N load. -Used in Super Tilt Reclining Mechanism. -Used in Low Rider, Reclining Mechanism. -Used in Low Rider, Both tilting cylinders (<20" seat width). -Used in MLT700, Tilting Mechanism. -Used in MLT700, Reclining Mechanism.		
TRP004 G600	600 N cylinder Long Sold with [CC009]. (SUSPA cylinder Trigger)	Long Suspa cylinder having a 600 N load. -Used in Low Rider, Both tilting cylinders (>20" seat width).		
TRP004 G150	150 N cylinder Long Sold with [CC009]. (SUSPA cylinder Trigger)	Long Suspa cylinder having a 150 N load. -Used in MLT700 footrest elevating system.		



TRP004 G200S	200 N cylinder Short Sold with [CC009]. (SUSPA cylinder Trigger)	Short Suspa cylinder having a 200 N load. -Used in Jr. Tilt Commode, Primary Tilt Cylinder. -Used in Tilt Commode (>18" seat width), Both Tilting Cylinders. -Used in NRG+ Recline (<22" seat width) Reclining Mechanism. -Used in Super Tilt, Secondary Tilt Cylinder (<24" seat width). -Used in MLTR footrest elevating system.	220/ each	
TRP004 G500S	500 N cylinder Short Sold with [CC009]. (SUSPA cylinder Trigger)	Short Suspa cylinder having a 500 N load. -Used in Super Tilt, Primary Tilt Cylinder (<24" seat width). -Used in MLTR, Tilting mechanism. -Used in MLTR, Reclining mechanism.		
TRP004 G400S	400 N cylinder Short Sold with [CC009]. (SUSPA cylinder Trigger)	Short Suspa cylinder having a 400 N load. -Used in NRG+ Recline (>22" seat width) Reclining Mechanism. -Used in Super Tilt, Secondary Tilt Cylinder (>24" seat width).		
TRP004 G800S	800 N cylinder Short Sold with [CC009]. (SUSPA cylinder Trigger)	Short Suspa cylinder having an 800 N load. -Used in Super Tilt, Primary Tilt Cylinder (>24" seat width).		

SCO	Suspa Cylinder Trigger	Trigger that activates the cylinder mechanism.	10/ each	
CC010	Faucher Damper Cylinder	Faucher Damper Cylinder. -Used in Tilt Commode (16"-18" seat width), Secondary Damper Cylinder.	180/ each	
TRP00 40	Faucher 200 N cylinder	200 N Faucher cylinder used in Articulating footrest. -Used in the Articulating footrest mechanism.	200/ each	
CC012	Faucher Trigger	Trigger for the Faucher 200N Cylinder	10/ each	
CC013	Faucher Head	Cylinder Head for the Faucher 200N Cylinder	10/ each	

MLPTACL	Electric Tilt Cylinder Power Tilt Dynamic Tilt	Motor activated Tilt Cylinder for the Power Tilt and the Dynamic Tilt. (Chair Serial number needs to be mentioned)	240/ each	
MLPTCBL	Battery Pack Power Cylinder Power Tilt Dynamic Tilt	Battery Pack for the Power Cylinder. (Used in the Power Tilt and the Dynamic Tilt)	40/ each	
STCC	Power Cylinder Charger Set Power Tilt Dynamic Tilt	Charger Set for the Power Cylinder. (Used in the Power Tilt and the Dynamic Tilt)	40/ each	
STM	Power Tilt Activation Switch Power Tilt Dynamic Tilt	Activation Switch for the Power Cylinder. (Used in the Power Tilt and the Dynamic Tilt)	60/ each	

CABG002	Gas Cylinder Cable Dual Head	<p>Cable for the Suspa gas cylinder having a Dual head.</p> <p>-Used in Super Tilt and Low Rider tilting mechanism (Seat width <20").</p> <p>-Used in MLT700 Tilt mechanism (Seat width <20").</p> <p>-Used in Tilt Commode, Tilting mechanism (Seat width 18"-20").</p>	80/ each	
CABG002 -W	Gas Cylinder Cable Dual Head (Long)	<p>Long Cable for the Suspa gas cylinder having a Dual head.</p> <p>-Used in Super Tilt and Low Rider tilting mechanism (Seat width >20").</p> <p>-Used in MLT700 Tilt mechanism (Seat width >20").</p> <p>-Used in Tilt Commode, Tilting mechanism (Seat width >20").</p>		

CABG001	Gas Cylinder Cable Single Head	<p>Standard Cable for the Suspa gas cylinder having a single head.</p> <p>-Used in NRG+ Recline Reclining mechanism (Seat width<20").</p> <p>-Used in Super Tilt Recline mechanism (Seat width<20").</p> <p>-Used in Low Rider Recline mechanism (Seat width<20").</p> <p>-Used in MLTR Tilting mechanism (Seat width<20").</p> <p>-Used in MLTR Elevating Footrest mechanism</p> <p>-Used in MLT700 Elevating Footrest mechanism</p> <p>-Used in Tilt Commode Tilting mechanism (Seat width 16"-18").</p>	60/ each	
CC052	Gas Cylinder Cable Single Head (Long)	<p>Long Cable for the Suspa gas cylinder having a single head.</p> <p>-Used in NRG+ Recline Reclining mechanism (Seat width<20").</p> <p>-Used in Super Tilt Recline mechanism (Seat width<20").</p> <p>-Used in Low Rider Recline mechanism (Seat width<20").</p> <p>-Used in MLTR Tilting mechanism (Seat width<20").</p>		

TRCH006	Cable Release Handle	Tilt and Recline release handle for all tilting and reclining wheelchairs.	30/ each	
PW010R	Platform Walker Brake Complete Assembly R	Right Side Brake assembly for the Platform Walker.	100/ set	
PW010L	Platform Walker Brake Complete Assembly L	Left Side Brake assembly for the Platform Walker.		
BC	Platform Walker Brake Cable and Lever	Brake Lever and Cable for the Platform Walker.	80/ each	