

Maple Leaf WHEELCHAIR






PARTS CATALOGUE Frames









6540 Tomken Road
Mississauga ON L5T 2E9
905.564.2250
905.564.5424
1.800.362.0873
www.mapleleafwheelchair.ca





Frames

Part No.	Part Name	Part Description	Cost	
SFA611	Access Standard Side Frame R	Standard Side frame for the Access Wheelchair for the Right-Hand Side. (Chair Serial Number needs to be mentioned)	440/ each	
SFA612	Access Standard Side Frame L	Standard Side frame for the Access Wheelchair for the Left-Hand Side. (Chair Serial Number needs to be mentioned)		
SFA613	Access Hemi-Side Frame R	Hemi Side frame for the Access Wheelchair for the Right-Hand Side. (Chair Serial Number needs to be mentioned)	440/ each	
SFA614	Access Hemi-Side Frame L	Hemi Side frame for the Access Wheelchair for the Left-Hand Side. (Chair Serial Number needs to be mentioned)		
SFSW607	Swift Standard Side Frame R	Standard Side frame for the Swift Wheelchair for the Right-Hand Side. (Chair Serial Number needs to be mentioned)	520/ each	
SFSW608	Swift Standard Side Frame L	Standard Side frame for the Swift Wheelchair for the Left-Hand Side. (Chair Serial Number needs to be mentioned)		




(All Images are for Visual reference please use the Part Description for determining the correct required part.)
Maple Leaf Wheelchair © 2021





SFSW609	Swift Hemi-Side Frame R	Hemi- Side frame for the Swift Wheelchair for the Right-Hand Side. (Chair Serial Number needs to be mentioned)	520/ each	
SFSW610	Swift Hemi-Standard Side Frame L	Hemi-Standard Side frame for the Swift Wheelchair for the Left-Hand Side. (Chair Serial Number needs to be mentioned)		
SFSW612	Swift Standard Angle Side Frame R	Standard Angle Adjustable Side frame for the Swift Wheelchair for the Right-Hand Side. (Chair Serial Number needs to be mentioned)	520/ each	
SFSW613	Swift Standard Angle Side Frame L	Standard Angle Adjustable Side frame for the Swift Wheelchair for the Left-Hand Side. (Chair Serial Number needs to be mentioned)		
SFSW614	Swift Hemi-Angle Adjustable Side Frame R	Hemi-Angle Adjustable Side frame for the Swift Wheelchair for the Right-Hand Side. (Chair Serial Number needs to be mentioned)	520/ each	
SFSW615	Swift Hemi-Angle Adjustable Side Frame L	Hemi-Angle Adjustable Side frame for the Swift Wheelchair for the Left-Hand Side. (Chair Serial Number needs to be mentioned)		




SFN601	NRG+ Standard Side Frame R	Standard Side frame for the NRG+ Wheelchair for the Right-Hand Side. (Chair Serial Number needs to be mentioned)	480/ each	
SFN602	NRG+ Standard Side Frame L	Standard Side frame for the NRG+ Wheelchair for the Left-Hand Side. (Chair Serial Number needs to be mentioned)		
SFN605	NRG+ Hemi-Side Frame R	Hemi- Side frame for the NRG+ Wheelchair for the Right-Hand Side. (Chair Serial Number needs to be mentioned)	480/ each	
SFN606	NRG+ Hemi-Side Frame L	Hemi- Side frame for the NRG+ Wheelchair for the Left-Hand Side. (Chair Serial Number needs to be mentioned)		
SFNR615	NRG+ Recliner Standard Side Frame R	Standard Side frame for the NRG+ Recliner Wheelchair for the Right-Hand Side. (Chair Serial Number needs to be mentioned)	520/ each	
SFNR616	NRG+ Recliner Standard Side Frame L	Standard Side frame for the NRG+ Recliner Wheelchair for the Left-Hand Side. (Chair Serial Number needs to be mentioned)		





SFNR617	NRG+ Recliner Hemi- Side Frame R	Hemi- Side frame for the NRG+ Recliner Wheelchair for the Right-Hand Side. (Chair Serial Number needs to be mentioned)	520/ each	
SFNR618	NRG+ Recliner Hemi- Side Frame L	Hemi- Side frame for the NRG+ Recliner Wheelchair for the Left-Hand Side. (Chair Serial Number needs to be mentioned)		
SFEZ620	EZ Ride, NRG+ Gold Standard Side Frame	Standard Side frame for the EZ Ride, NRG+ Gold Wheelchair. (Chair Serial Number needs to be mentioned)	560/ each	
SFEZ621	EZ Ride, NRG+ Gold Hemi- Side Frame	Hemi- Side frame for the EZ Ride, NRG+ Gold Wheelchair. (Chair Serial Number needs to be mentioned)	560/ each	
SF621	(NRG+, Swift) Adjustable Height Standard Back Post	Height Adjustable Standard Back Post for the NRG+ and Swift. (Chair Serial Number needs to be mentioned)	140/ each	
SF622	(NRG+, Swift) Adjustable Height Tall Back Post	Height Adjustable Tall Back Post for the NRG+ and Swift. (Chair Serial Number needs to be mentioned)	160/ each	



(All Images are for Visual reference please use the Part Description for determining the correct required part.)
Maple Leaf Wheelchair © 2021





SF623	(NRG+ Gold, EZ Ride) Adjustable Height Standard Back Post	Height Adjustable Standard Back Post for the NRG+ Gold and EZ Ride. (Chair Serial Number needs to be mentioned)	140/ each	
SF624	(NRG+ Gold, EZ Ride) Adjustable Height Tall Back Post	Height Adjustable Tall Back Post for the NRG+ Gold and EZ Ride. (Chair Serial Number needs to be mentioned)	160/ each	
MLREC	NRG+ Recliner Back Post	Standard Back Post for the NRG+ Recliner. (Chair Serial Number needs to be mentioned)	180/ each	
SF5SN	(NRG+, Swift) Telescopic Back Post Sleeve	Telescopic Back Post Sleeve for the NRG+ and Swift. (Chair Serial Number needs to be mentioned)	75/ each	
SF5EG	(NRG+ Gold, EZ Ride) Telescopic Back Post Sleeve	Telescopic Back Post Sleeve for the NRG+ Gold and EZ Ride. (Chair Serial Number needs to be mentioned)	75/ each	

UPHA519	NRG+ Recliner Headrest Tube	Headrest Tube for the NRG+ Recliner Mode. (Chair Serial Number needs to be mentioned)	55/ each	
FRM209	NRG+ Recliner Post Bracket Set	NRG+ Recliner foldable back post bracket set. (Chair Serial Number needs to be mentioned) Sold with 2x Screw [HDW033], 2x Thin Nut [HDW009] 2x Custom Mounting Bushing [BB019]	120/ set	
SFST701	Power Super Tilt Lower Base	Lower Base of the Super Tilt Wheelchair. (Chair Serial Number needs to be mentioned)	480/ each	
SFST701-ST	Super Tilt Lower Base for Power Tilt	Lower Base of the Super Tilt Wheelchair having Power Tilt option. (Chair Serial Number needs to be mentioned)	520/ each	
SFST702	Super Tilt Upper Frame	Upper Frame of the Super Tilt Frame. (Chair Serial Number needs to be mentioned)	360/ each	
SFST702 ST	Super Tilt Upper Frame for Power Tilt	Upper Frame of the Super Tilt Frame having Power Tilt option. (Chair Serial Number needs to be mentioned)	400/ each	





SFST701-D	Low Rider Lower Base	Lower Base of the Low Rider Wheelchair. (Chair Serial Number needs to be mentioned)	480/ each	
FRM022P	Low Rider Lower Base for Power Tilt	Lower Base of the Low Rider Wheelchair having Power Tilt option. (Chair Serial Number needs to be mentioned)	520/ each	
SFST702-D	Low Rider Upper Frame	Upper Frame of the Low Rider Frame. (Chair Serial Number needs to be mentioned)	360/ each	
FRM023P	Low Rider Upper Frame for Power Tilt	Upper Frame of the Low Rider Frame having Power Tilt option. (Chair Serial Number needs to be mentioned)	400/ each	
SFST710	Super Tilt and Low Rider Back Post	Standard Back Post for the Super Tilt and Low Rider. (Chair Serial Number needs to be mentioned)	140/ each	
SFST711	Super Tilt and Low Rider Straight Back Post	Straight Back Post for the Super Tilt and Low Rider. (Chair Serial Number needs to be mentioned)	140/ each	
SFST710 EM	Super Tilt and Low Rider Endomorphic Back Post	Endomorphic wider Back Post for the Super Tilt and Low Rider. (Chair Serial Number needs to be mentioned)	140/ each	





FRM027	Super Tilt and Low Rider Recline Back post	<p>Recline back post with welded upper crosstie used in Super Tilt and Low Rider wheelchairs.</p> <p>(Chair Serial Number needs to be mentioned)</p>	280/ each	
FRM123	Low Rider Frame Support Bracket	Bracket that sits on Lower base that supports the Upper frame for the Low Rider wheelchair.	50/ each	
FRM121	Super Tilt Frame Support Bracket	Bracket that sits on Lower base that supports the Upper frame for the Super Tilt wheelchair.	50/ each	
SFST716	Super Tilt and Low Rider Recline Lower Cross Tie	<p>Cross Tie for the Super Tilt and Low Rider Recline Wheelchair.</p> <p>(Chair Serial Number needs to be mentioned)</p>	140/ piece	

SFST712	7/8" Stroller Handler Super Tilt Low Rider	Stroller handle for the Super Tilt and Low Rider wheelchairs. Sold with 2x [PPG013]. (Wheelchair Serial number needs to be mentioned)	140/ piece	
SFSW712	1" Stroller Handler (Steel) Swift NRG+ NRG+ Recline	Stroller handle for the Swift and NRG+ wheelchairs. Sold with 2x [PPG013]. (Wheelchair Serial number needs to be mentioned)	140/ piece	
RSH	1" Stroller Handler (Aluminum) NRG+ Gold EZ Ride	Stroller handle for the NRG+ Gold and EZ Ride wheelchairs. Sold with 2x [PPG013]. (Wheelchair Serial number needs to be mentioned)	140/ piece	
ASHJ	Adjustable Stroller Handle Joint for Super Tilt and Low Rider	Adjustable Stroller Handle Joint used to mount Stroller Handle [FRM012] to Back Post on Super Tilt and Low Rider wheelchair models. Sold with the Pair of 2x Screw [HDW092] [#10-24 x 1" Long Phillips-Square Machine Screw] 2x Nut [HDW069] [#10-24 x 13mm Joint Connection Nut (UNC)]	70/ each	




STFST712 RSHHC	7/8" Foam Grip Spreader Bar Super Tilt Low Rider	7/8" Foam grip for the spreader bar used on the Super Tilt and the Low Rider wheelchair.	10/ each	
FG	1" Foam Stroller Handle Grip	1" Foam grip for the stroller handles for the Swift, NRG+, NRG+ Recline, NRG+ Gold and the EZ Ride wheelchair models.		
UPHA514 -4	3/4" Foam Back Post Grip. NRG+ Gold EZ Ride	3/4" Foam grip for the NRG+ Gold and EZ Ride Back Post.	10/ each	
UPHA514	7/8" Foam Back Post Grip. Swift NRG+ NRG+ Recline	7/8" Foam grip for Swift, NRG+ and NRG+ Recline Back Post.		
SFML801	MLTR 500 Lower Base	Lower Base for the MLTR 500 Wheelchair. (Chair Serial Number needs to be mentioned)	440/ each	
SFML803	MLTR 600 Lower Base	Lower Base for the MLTR 600 Wheelchair. (Chair Serial Number needs to be mentioned)	440/ each	





SFML804	MLTR Upper Frame	Upper Frame for the MLTR 500/600 Wheelchair. (Chair Serial Number needs to be mentioned)	520/ each	
SFML809	MLTR Lower Back Frame	Lower Back Frame for the MLTR 500/600 Wheelchair. (Chair Serial Number needs to be mentioned)	280/ each	
SFML810	MLTR Upper Back Frame	Upper Back Frame for the MLTR 500/600 Wheelchair. (Chair Serial Number needs to be mentioned)	320/ each	
UPHML311	MLTR Head Wing mount	Head Wing mount for the MLTR.	80/ each	






SFML805	MLT 700 Lower Base	Lower Base for the MLT 700 A/B Wheelchair. (Chair Serial Number needs to be mentioned)	480/ each	
SFML806	MLT 700 Upper Frame	Upper Frame for the MLT 700 A/B Wheelchair. (Chair Serial Number needs to be mentioned)	320/ each	
SFML702	MLT 700B Frame Height Adjuster	Tubular frame Height Adjuster used in MLT 700B. Used to level front caster wheels to the back. Sold with 2x Bolt [B2H12] 2x Nut [HDW009] 4x Washer [HDW070]	70/ each	
UPHT226 TBO	MLT 700 Head Wing mount	Head Wing mount for the MLT700, upholstery sits on mount.	80/ each	






SHML807	MLT 700 Back Frame	<p>Back Frame of the MLT 700 A/B Wheelchair.</p> <p>(Chair Serial Number needs to be mentioned)</p>	280/ each	
UPHT226 TBO	MLT 700 Head Wing mount	Head Wing mount for the MLT700, upholstery sits on mount.	80/ each	
HBAL	<p>Prism-Block Clamp Assembly Type 1</p> <p>(Square-Circle)</p>	<p>Prism-Block shaped Clamp assembly that attaches head wing to frame for MLT700. (Square-Circle)</p> <p>Sold with [HKH]</p>	55/ each	
HKH	Prism-Block Clamp Knob	Knob used to tighten both Prism Block pieces.	15/ each	

(All Images are for Visual reference please use the Part Description for determining the correct required part.)
Maple Leaf Wheelchair © 2021







HDC001	HD Commode Frame	Complete Frame Assembly for the HD Commode COM1. (Chair Serial Number needs to be mentioned)	1000/ each	
HDC002	HD Commode Frame	HD Commode with Heavy Duty COM2.	1100/ each	
HDC003	HD Commode Frame	HD Commode with Heavy Duty COM3.	1200/ each	
TSC001	Tilt Commode Lower Base	Lower Base for the Tilt Commode. (Chair Serial Number needs to be mentioned)	520/ each	
TSC002	Tilt Commode Upper Frame	Upper Frame for the Tilt Commode. (Chair Serial Number needs to be mentioned)	660/ each	

JTSC001	Jr. Tilt Commode Lower Base	Lower Base for the Jr. Tilt Commode. (Chair Serial Number needs to be mentioned)	520/ each	
JTSC002	Jr. Tilt Commode Upper Frame	Upper Frame for the Jr. Tilt Commode. (Chair Serial Number needs to be mentioned)	660/ each	
PPG005	Refresh Commode Quick Release Pin	Quick Release Pin attaches the Refresh Commode Bracket (FRM318) to the Refresh Commode Base (FRM316). [5/16" x 1 1/8"]	10/ each	
FRM901	Refresh Commode Base	Lower Base for the Refresh Commode. (Chair Serial Number needs to be mentioned)	550/ each	


FRM902	Refresh Commode Frame	Upper Frame for the Refresh Commode. (Chair Serial Number needs to be mentioned)	650/ each	
FRM903	Refresh Commode Bracket	Refresh Commode Bracket that attaches base to the upper frame.	40/ each	
FRM904	Refresh Commode Cylinder Bracket	Refresh Commode Cylinder Bracket houses cylinder that connects frame and base.	40/ each	
HBAL	Prism-Block Clamp Assembly Type 1 (Square-Circle)	Prism-Block shaped Clamp assembly that attaches Lateral to frame for Jr. Tilt Commode, Tilt Commode, and the Refresh Commode. (Square-Circle) Sold with [PPG007]	55/ each	
HKH	Prism-Block Clamp Knob	Knob used to tighten both Prism Block pieces.	15/ each	

HWLB-R	Commode Lateral R	Right Side Lateral mounting bracket used for Tilt, Jr. Tilt and Refresh Commode.	100/ each	
HWLB-L	Commode Lateral L	Left Side Lateral mounting bracket used for Tilt, Jr. Tilt and Refresh Commode.	100/ each	
FRM601R	Escape Commode Side Frame R	Right Side frame for the Escape Commode. (Chair Serial Number needs to be mentioned)	550/ each	
FRM602L	Escape Commode Side Frame L	Left Side frame for the Escape Commode. (Chair Serial Number needs to be mentioned)		
PW009	Adult upright frame	Adult upright frame for the Platform Walker.	180/ each	
PW010	Junior upright frame	Junior upright frame for the Platform Walker.	170/ each	
PW007	Adult base	Adult base frame for the Platform Walker.	220/ each	
PW008	Junior base	Junior base frame for the Platform Walker.	210/ each	

(All Images are for Visual reference please use the Part Description for determining the correct required part.)
Maple Leaf Wheelchair © 2021

PW005	Adult upright brace	Adult upright brace for the Platform Walker.	200/ each	
PW006	Junior upright brace	Junior upright brace for the Platform Walker.	190/ each	
PW004	Height adjustment bracket	Height adjustment bracket for the platform cushion. Sold as a pair with 8x Mounting screws. (HDW036)	60/ each	
PW003	Arthritic Handle with grips	Handle for the Platform Walker.	120/ each	
SFML815	Pivot Bushing	Delrin bushing for the MLTR pivot and Tilt and Jr. Tilt Commode. (Hinge where wheelchair tilts)	10/ each	
TSC007	NRG+ Recline Post Bracket Bushing	Bushing on the NRG+ Recline Post that pivots the folding bracket and the frame.	10/ each	

BB020	Refresh Commode Cylinder Bushing	Cylinder Bushing for the Refresh Commode.	10/ each	
STCLMP 002	7/8" Cable Clamp	Cable Clamp for the 7/8" Tubing.	5/ each	
STCLMP 001	1" Cable Clamp	Cable Clamp for the 1" Tubing.	5/ each	
HWF003	1" Standard Plug	Plug used for all standard 1" tubing.	1/ each	
HWF002	7/8" Standard Plug	Plug used for all standard 7/8" tubing.	1/ each	
HWF001	3/4" Standard Plug	Plug used for all standard 3/4" tubing.	1/ each	

HWF010	5/8" Standard Plug	Plug used for all standard 5/8" tubing. (Used in Head Wing Lateral for commode.)	1/ each	
HWF007	3/4" Square Universal Plug	Square Plug for all chairs with 3/4" square piping. (Used in wheelchair frames)	1/ each	
HWF008	7/8" Square Universal Plug	Square Plug for all chairs with 7/8" square piping. (Used in wheelchair frames)	1/ each	
HWF009	1" Square Universal Plug	Square Plug for all chairs with 1" square piping. (Used in wheelchair frames)	1/ each	
PPG033	Super Tilt Oval Plug	Oval Plug for the Super Tilt Wheelchair.	1/ each	